



SYMMETRIKIT
AVANTI®

User Handbook



SYMMETRIKIT
POSTURAL CARE

HH a division of The Helping Hand Company

-
- *Your new chair* 2
 - *Set up checklist* 4
 - *Getting to know your chair* 5-6
 - *Looking after your chair* 6
 - *Guarantee and warranty* 7
 - *How to get in touch with us* 7
 - *Frequently asked questions* 8



YOUR NEW CHAIR

- Your AVANTI chair is adjusted to your personal settings, keeping you comfortable and supported all day long. This ensures peace of mind every day to be able to get up and about. No slumping or perching on the edge of the seat as your feet are safely anchored to the floor.
- Sit right every time – as the chair is adjusted to fit you, each time you sit down you will be positioned to the back of the chair with no need to shuffle back or reposition to get comfortable. Maintenance of good posture and back support when seated.
- Multi adjustable - the back, seat and leg support angle can be individually adjusted to create your own snooze zone.

- ✔ Adjustable postural support
- ✔ Medium / high pressure relief as standard
- ✔ Mobile and semi ambulant
- ✔ Maintaining independence
- ✔ Tilt in space, recline and leg rest elevation functions as standard
- ✔ Up to 150kg/24 stone
- ✔ Moving room to room
- ✔ Single point easy brake system

Your chair - a perfect fit for you

Back rest with optional built in lateral support to help maintain good posture

Single point easy brake system

Adjustable recline and tilt in space.

Medium / high pressure management as standard

Height adjustable leg rest, raise up to 45°



SET UP CHECKLIST

We've included a helpful checklist for delivery of your chair

Home access

- Access to main entrance (ramp/steps) clear of obstructions/trip hazards for delivery vehicle
(Home delivery via two person team as standard)
- External doorway wide enough, minimum doorway width 730mm
- Confirmed access if delivery involves use of lift or carry up in internal stairs

Room setting

- Internal door access width, minimum 730mm
- Adequate space around chair to allow back to fully recline
- Each function of chair checked and working
- Demonstration on how to use chair
- Test function and fit of chair

and relax...

Note: we have completed this for you during our final quality check

Serial Number:

Date of Manufacture:

The serial number is located on the rear frame of the chair.
It is important to reference the serial number during all communication with your local supplier/The Helping Hand Company.

GETTING TO KNOW YOUR CHAIR

Tilt in space



- Pull the right hand lever on the back of the seat to select **TILT** positions of the chair

Back recline



- Pull the left hand lever on the back of the seat to select back **RECLINE** positions of the chair

Easy brake system



- Pull the bar when you want to move the chair - let it go and both rear castors immediately lock keeping it stable when standing/sitting



 <p>Foot plate</p>	<ul style="list-style-type: none">• To lower the foot plate, release the lever underneath as shown in the left photo	 <p>Leg rest</p>	<ul style="list-style-type: none">• To adjust the angle, use the lever located on the top right hand side of the leg rest
------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------

LOOKING AFTER YOUR CHAIR

- Dust and dirt can increase the wear and tear on the fabric, we recommend you brush/vacuum the chair regularly
- If the upholstery is exposed to direct sunlight for long periods of time it will likely fade, where possible use blinds/curtains for protection
- Cuero upholstery should be wiped down using a damp, clean cloth with a solution of mild detergent and lukewarm water
- The frame of the chair can also be cleaned with a damp cloth and mild detergent if required.

GUARANTEE AND WARRANTY

Mechanical parts - two years

Upholstery - one year

- Normal wear and tear, soiling of covers and fabric shrinkage are excluded. We reserve the right to withdraw colours or materials from the range but will always endeavour to find the closest match available
- This chair is not intended for use outdoors OR within a vehicle
- Great care is taken in the design and manufacture of your chair and to ensure it meets all safety standards - do not allow anyone to kneel, stand on the leg rest, arms or back, as incorrect use of this product can cause accidents



HOW TO GET IN TOUCH WITH US



Tel: +44 (0) 1531 635388
Email: sales@Symmetrikit.com
Online: www.Symmetrikit.com

FREQUENTLY ASKED QUESTIONS

- **What if I have a problem with my chair?**

In the highly unlikely event of a problem, please refer to page 5 - 7

- **Can I move from room to room?**

*The chair can be moved with someone sitting in it as it is designed to be moved around the house.
It should not be used outdoors*

- **How do I keep my chair in tip top condition?**

*We recommend an annual inspection to check for wear and tear, and test mechanical integrity.
For further details please contact our Service team on 01531 634468 - remember to have your serial number handy.*

- **How do I go about obtaining new upholstery?**

We reserve the right to withdraw colours or materials from the range but will always endeavour to find the closest match available. Please have the serial number of the chair to hand when you call.

- **How will you use my personal details?**

*Your details are stored on a secure database and are protected by the Data Protection Act 1984.
We only store your details to cover your guarantees etc, we never pass your details to third parties.
Credit or debit card details are never stored and online card payment details are processed by a regulated merchant service to protect you.*



The Helping Hand Company
Bromyard Road, Ledbury
Herefordshire HR8 1NS

Tel: 01531 635678
www.AVANTChair.co.uk
sales@HelpingHand.co.uk

Tel: 01531 635388
www.Symmetrikit.com
sales@Symmetrikit.com



SKLX9028_2



HH a division of The Helping Hand Company