

ROA011R - Anti Roll-Back Lock Set

These can be engaged or disengaged using the spring pull pin. When engaged the one-way knurled bearing assembly will rest on the wheel and allow the unit to roll forward but not backwards when in use.



ROA007LR - Large Headrest

Similar in function to the standard headrest except the overall width is increased to 230mm.



ROA007SR - Special Headrest

A special headrest can be made to order to suit the end user.



ROA012R - Central Positioning Bar with Abduction Pad

Mounted off the profile extrusion in place of the standard central positioning bar. It can be used to help prevent leg scissoring.





Moorings Mediquip 51 Slaght Road Ballymena Co. Antrim United Kingdom BT42 2JH

ROA008 - Swing away shoulder positioners

Mounted off a horizontal bar attached to the front profile extrusion. These are adjustable in height, width and angle to deal with different size variations. 2 sizes of pad available.



ROA009 - Handholds (Straight)

Mounted from the base frame and provide an armrest and hand grip handle for the client. The framework is adjustable in height and length to grow with the user, the grip is also angle adjustable.



ROA009B - Handhold (Bumper Style)

Mounted off the handhold uprights to act as a simple rail for the user to hold during use, with Vertical and horizontal adjustment. The addition of a plastic bracket to the grab rail can also be used to hold an iPad.



