



Specifications

Frame size:	Short	Tall
Width	???mm	???mm
Depth	???mm	???mm
Leg length	???mm	???mm
Angle back	???°	???°
Angle legrest	???°	???°
Angle footplate	???°	???°
Tilt	???°	???°
Seat height	???mm	???mm
Overall length tilted & leg up		
Manual tilt, recline & leg angle operation		
Weight	???Kg	???Kg



©HHCo ©Symmetrikit 2018 - 090818



Simple supportive seating



Easy Adjustment

Printed on recycled paper with vegetable inks



The Helping Hand Company
Bromyard Road, Ledbury
Herefordshire HR8 1NS

Tel: 01531 635388
www.Symmetrikit.com
sales@Symmetrikit.com



HH a division of The Helping Hand Company

SYMMETRIKIT AVANTI®

Simple supportive seating with built in pressure care*



Quick and easy adjustment makes it easy to get the seat depth and width right. It's a great start to reducing point pressures as well as stabilising the seating base.

Setting the correct leg length and arm height helps promote stability and further spreads the load.



Tilt in space

Forward for function or access
Back to offload the bony prominences and to relax.

An independently adjustable back

Allows you to get the final fit right.

***Lowzone™ Contour cushion as standard**

Our contoured Lowzone™ visco-elastic cushions are regarded as providing "Medium High" pressure relief. In our own trials using our FSA pressure monitors, these regularly exceed non contoured gel, water and standard visco-elastic cushions in pressure reduction. They are inherently more stable than fluid based systems.

Options

2 Frame heights:

Short for ambulant users or Tall for taller, semi ambulant or sling lifted users.

Laterals for extra support

Head cushion

Footplate for non ambulant users

Tray

Single point braking

Tool free adjustment

can be switched for hex keys where no adjustment is required.

